


patient's ID : demo-ERG-moncolor
file number :
birth date : 23/10/1987


Rx :
exam. date : 11/09/2009 08:33
exam :

STANDARD ISCEV ERG EXAM WITH ERG JET SENSOR


Scotopic 0.01 ERG 19mn 28s Val= 3 Rej= 0
BI stimulated




Scotopic 3.0 ERG 24mn 43s Val= 5 Rej= 0
BI stimulated




Scotopic 10 ERG 25mn 52s Val= 6 Rej= 0
BI stimulated



Scotopic 3.0 OPs 27mn 25s Val= 5 Rej= 0
BI stimulated



Photopic 3.0 ERG 37mn 19s Val= 6 Rej= 0
BI stimulated



Photopic 3.0 flicker 38mn 10s Val= 17 Rej= 0
BI stimulated

